

# Comparative performance of plasma A $\beta$ 42/A $\beta$ 40 and p-tau181 for the detection of early brain amyloid deposition in individuals with subjective cognitive decline

María Pascual-Lucas<sup>1</sup>, José Antonio Allué<sup>1</sup>, Leticia Sarasa<sup>1</sup>, Noelia Fandos<sup>1</sup>, Sergio Castillo<sup>1</sup>, Jose Terencio<sup>1</sup>, Manuel Sarasa<sup>1†</sup>, Juan Pablo Tartari<sup>2</sup>, Ángela Sanabria<sup>2,3</sup>, Lluís Tárraga<sup>2,3</sup>, Agustín Ruiz<sup>2,3</sup>, Marta Marquie<sup>2,3</sup>, Mercè Boada<sup>2,3</sup>

<sup>1</sup>Araclon Biotech-Grifols - Zaragoza (Spain), <sup>2</sup>Ace Alzheimer Center Barcelona - Universitat Internacional de Catalunya - Barcelona (Spain), <sup>3</sup>CIBERNED, National Institute of Health Carlos III - Madrid (Spain), <sup>†</sup> Deceased May 27, 2020



## BACKGROUND

In the last years, **blood-based biomarkers** have shown high accuracy for the identification of early alterations of Alzheimer's disease (AD). Among them, plasma A $\beta$ 42/A $\beta$ 40 and p-tau181 have been proved to be reliable biomarkers of brain amyloid deposition. However, the quantification of these molecules in plasma, mainly A $\beta$ 42, presents analytical difficulties and **high reliability analytical assays** are needed in order to avoid biased conclusions about the comparative performance of both biomarkers. In this study, we aim to avoid these uncertainties by the introduction of a high sensitivity assay based on HPLC-MS for the quantification of A $\beta$  peptides in plasma.

## OBJECTIVE

To compare the ability of **plasma A $\beta$ 42/A $\beta$ 40 ratio** measured with a mass spectrometry-based assay and **plasma p-tau181** measured with a high-sensitivity technology to detect **early brain amyloid deposition** in individuals at risk of AD.

## METHODS

152 subjects with subjective cognitive decline (SCD) from the FACEHBI cohort were included in the present study, of which 16% were A $\beta$ -PET (+). Plasma A $\beta$ 40 and A $\beta$ 42 were quantified with a high sensitivity antibody-free mass spectrometry-based assay (ABtest-MS, Araclon Biotech). Plasma p-tau181 was measured with Simoa® pTau-181 V2 Advantage Kit (Quanterix). The ability of plasma biomarkers alone, combined or after the addition of demographic covariates, to detect A $\beta$ -PET positivity was assessed by logistic regression and ROC curve analysis.


## RESULTS

### 1. Characteristics of the study population

**Table 1. Characteristics of the study population.** Data are median values (interquartile range) or number of cases (%). Differences between A $\beta$ -PET (-) and A $\beta$ -PET (+) groups were tested using Mann-Whitney and Chi-square tests, as appropriate.

|                           | All                 | A $\beta$ -PET (-)  | A $\beta$ -PET (+)  | P value |
|---------------------------|---------------------|---------------------|---------------------|---------|
| Participants              | 152                 | 125 (82%)           | 27 (18%)            |         |
| Age, years                | 66.0 (60.0-70.0)    | 64 (60.0-69.0)      | 70 (67.0-72.0)      | .001    |
| Female                    | 95 (63%)            | 82 (66%)            | 13 (48%)            | .139    |
| APOE $\epsilon$ 4         |                     |                     |                     | <.0001  |
| 1 allele                  | 39 (26%)            | 23 (18%)            | 16 (59%)            |         |
| 2 alleles                 | 3 (2%)              | 2 (2%)              | 1 (4%)              |         |
| MMSE, score               | 30 (29-30)          | 29 (29-30)          | 30 (29-30)          | .361    |
| FBB-PET, CL               | -1.43 (-6.66-6.57)  | -3.6 (-7.8-1.5)     | 32.8 (21.2-60.6)    | <.0001  |
| A $\beta$ 42/A $\beta$ 40 | 0.258 (0.236-0.276) | 0.261 (0.244-0.278) | 0.216 (0.199-0.236) | <.0001  |
| p-tau181, pg/ml           | 1.81 (1.38-2.40)    | 1.71 (1.32-2.07)    | 2.67 (2.28-3.19)    | <.0001  |

Plasma A $\beta$ 42/A $\beta$ 40 levels were significantly reduced in the A $\beta$ -PET (+) group compared to A $\beta$ -PET (-) subjects, whereas the opposite was observed for plasma p-tau181 ( $P$ <.0001 in both cases, Mann-Whitney test) (Table 1 and Figure 1).



**Figure 1. Box and whiskers plots of plasma A $\beta$ 42/A $\beta$ 40 and p-tau181 between A $\beta$ -PET groups.** \*\*\*  $P$  < .0001, Mann-Whitney test.


### 2. Discriminative ability of A $\beta$ -PET status

Plasma A $\beta$ 42/A $\beta$ 40 ratio and p-tau181 identified A $\beta$ -PET status with an AUC of 0.86 (95% CI 0.78–0.94) and 0.83 (95% CI 0.76–0.90), respectively (Figure 2). At the maximum Youden index, both plasma biomarkers presented a sensitivity of 81.5%, whereas A $\beta$ 42/A $\beta$ 40 ratio presented slightly superior specificity (84.0%) and overall accuracy (83.6%) than p-tau181 (80.8% and 80.9%, respectively) (Table 2).

The combination of plasma biomarkers yielded an AUC of 0.88 (95% CI 0.82–0.95). The inclusion of demographic covariates (age and APOE  $\epsilon$ 4 number of alleles) increased the AUC up to 0.90 (95% CI 0.85–0.96) (full model). Both A $\beta$ 42/A $\beta$ 40 and p-tau181 contributed significantly to this model ( $P$ <.0001 and  $P$ =.04, respectively). The full model significantly outperformed p-tau181 alone ( $\Delta$ AUC=0.07,  $P$ =.04), but did not differ from A $\beta$ 42/A $\beta$ 40 ( $\Delta$ AUC=0.04,  $P$ =.09) (Table 2).

Goodness of fit was assessed using AIC. Plasma A $\beta$ 42/A $\beta$ 40 model presented lower AIC (101.3) than the model composed of plasma p-tau181 (AIC=115.7).

The combination of both biomarkers yielded AIC=94.7, and the inclusion of demographic covariates did not further improve the fit in this case (AIC=94.4).



**Figure 2. ROC curves for discriminating A $\beta$ -PET status.**

**Table 2. Performance of predictive models to identify A $\beta$ -PET status.** \*  $P$  values correspond to the comparison of AUC versus the full model using DeLong test.

| Biomarker 1               | Biomarker 2 | Covariates             | AUC (95% CI)     | $P$ value*  | Sen (%) | Sp (%) | PPV (%) | NPV (%) | Acc (%) |
|---------------------------|-------------|------------------------|------------------|-------------|---------|--------|---------|---------|---------|
| A $\beta$ 42/A $\beta$ 40 | -           | -                      | 0.86 (0.78-0.94) | ns          | 81.5    | 84.0   | 52.4    | 95.5    | 83.6    |
| p-tau181                  | -           | -                      | 0.83 (0.76-0.90) | <b>0.04</b> | 81.5    | 80.8   | 47.8    | 95.3    | 80.9    |
| A $\beta$ 42/A $\beta$ 40 | p-tau181    | -                      | 0.88 (0.82-0.95) | ns          | 85.2    | 86.4   | 57.5    | 96.4    | 86.2    |
| A $\beta$ 42/A $\beta$ 40 | p-tau181    | Age, APOE $\epsilon$ 4 | 0.90 (0.85-0.96) | -           | 92.6    | 75.2   | 44.6    | 97.9    | 78.3    |

## CONCLUSIONS

The use of a high reliability method for the quantification of plasma A $\beta$  based on HPLC-MS (ABtest-MS) suggests that **plasma A $\beta$ 42/A $\beta$ 40 ratio could be a more accurate biomarker of early brain amyloid deposition** than p-tau181 in the first stages of AD.